

Corticosteroid 11-beta-dehydrogenase isozyme 2

Catalog No: #AP82310



Package Size: #AP82310-1 50ug #AP82310-2 100ug #AP82310-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Corticosteroid 11-beta-dehydrogenase isozyme 2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | 11-beta-hydroxysteroid dehydrogenase type 2,11-DH2,11-beta-hydroxysteroid dehydrogenase type II,11-HSD type II,NAD-dependent 11-beta-hydroxysteroid dehydrogenase,11-beta-HSD,Short chain dehydrogenase/reductase family 9C member 3,HSD11K,SDR9C3 |
| Accession No. | Uniprot ID:P80365 |
| Uniprot | P80365 |
| GeneID | 3291; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSD11B2

Note: This product is for in vitro research use only