

Catenin beta-1

Catalog No: #AP82313

Package Size: #AP82313-1 50ug #AP82313-2 100ug #AP82313-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Catenin beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Beta-catenin,Catnb
Accession No.	Uniprot ID:Q9WU82
Uniprot	Q9WU82
GeneID	84353;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ctnnb1

Note: This product is for in vitro research use only