

Beta-defensin 3

Catalog No: #AP82314



Package Size: #AP82314-1 50ug #AP82314-2 100ug #AP82314-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Beta-defensin 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	BNBD-3,BNDB-3,BNDB3
Accession No.	Uniprot ID:P46161
Uniprot	P46161
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DEFB3

Note: This product is for in vitro research use only