

High affinity copper uptake protein 1

Catalog No: #AP82315



Package Size: #AP82315-1 50ug #AP82315-2 100ug #AP82315-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	High affinity copper uptake protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Copper transporter 1,rCTR1,Liver regeneration-related protein LRRGT00200,Solute carrier family 31 member 1,Copt1,Ctr1
Accession No.	Uniprot ID:Q9JK41
Uniprot	Q9JK41
GeneID	171135;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc31a1

Note: This product is for in vitro research use only