

Protein CYR61

Catalog No: #AP82323



Package Size: #AP82323-1 50ug #AP82323-2 100ug #AP82323-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein CYR61
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CCN family member 1,Cysteine-rich angiogenic inducer 61,Insulin-like growth factor-binding protein 10,IBP-10,Protein GIG1,CCN1,GIG1,IGFBP10
Accession No.	Uniprot ID:O00622
Uniprot	O00622
GeneID	3491;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CYR61

Note: This product is for in vitro research use only