

## DNA (cytosine-5)-methyltransferase 1

Catalog No: #AP82324



Package Size: #AP82324-1 50ug #AP82324-2 100ug #AP82324-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DNA (cytosine-5)-methyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dnmt1,DNA methyltransferase Mmul,DNA MTase Mmul,MCMT,Dnmt,Met1,Uim
Accession No.	Uniprot ID:P13864
Uniprot	P13864
GeneID	13433;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Dnmt1

Note: This product is for in vitro research use only