

Deoxyribonuclease gamma

Catalog No: #AP82325



Package Size: #AP82325-1 50ug #AP82325-2 100ug #AP82325-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Deoxyribonuclease gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DNase gamma, DNaseY, Deoxyribonuclease I-like 3, DNase I-like 3
Accession No.	Uniprot ID:O89107
Uniprot	O89107
GeneID	116687;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Dnase113

Note: This product is for in vitro research use only