

Dihydropyrimidine dehydrogenase [NADP(+)]

Catalog No: #AP82326



Package Size: #AP82326-1 50ug #AP82326-2 100ug #AP82326-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dihydropyrimidine dehydrogenase [NADP(+)]
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DHPDHase,Dihydrothymine dehydrogenase,Dihydrouracil dehydrogenase
Accession No.	Uniprot ID:Q8CHR6
Uniprot	Q8CHR6
GeneID	99586;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dpyd

Note: This product is for in vitro research use only