

C-X-C chemokine receptor type 2

Catalog No: #AP82329



Package Size: #AP82329-1 50ug #AP82329-2 100ug #AP82329-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C chemokine receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CXC-R2,High affinity interleukin-8 receptor B,IL-8R B,CD182,IL8RB
Accession No.	Uniprot ID:Q28003
Uniprot	Q28003
GeneID	281863;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CXCR2

Note: This product is for in vitro research use only