

PAK4/5/6 (Phospho-Ser474) Antibody

Catalog No: #12417



Package Size: #12417-1 50ul #12417-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	PAK4/5/6 (Phospho-Ser474) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PAK4/5/6 (Phospho-Ser474) Antibody detects endogenous levels of PAK4/5/6 only when phosphorylated at Ser474
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PAK4/5/6 (Phospho-Ser474)
Target Name	PAK4/5/6
Modification	Phospho
Other Names	PAK4, KIAA1142, p21(CDKN1A)-activated kinase 4, PAK-4, p21-activated kinase 4
Accession No.	Swiss-Prot#: O96013NCBI Gene ID: 10298
Uniprot	O96013
GeneID	10298;
Target Species	human
Calculated MW	64kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:2000

Images

Western Blot analysis of lysates of HepG2 cell UV treated and HepG2 cell, using primary antibody at 1:1000 dilution.



Note: This product is for in vitro research use only