

C-X-C motif chemokine 16

Catalog No: #AP82330

Package Size: #AP82330-1 50ug #AP82330-2 100ug #AP82330-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Small-inducible cytokine B16,Transmembrane chemokine CXCL16
Accession No.	Uniprot ID:Q6AXU5
Uniprot	Q6AXU5
GeneID	497942;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cxcl16

Note: This product is for in vitro research use only