

Dual specificity protein phosphatase 2

Catalog No: #AP82333



Package Size: #AP82333-1 50ug #AP82333-2 100ug #AP82333-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dual specificity protein phosphatase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dual specificity protein phosphatase PAC-1,Pac-1,Pac1
Accession No.	Uniprot ID:Q05922
Uniprot	Q05922
GeneID	13537;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dusp2

Note: This product is for in vitro research use only