

Chymotrypsinogen 2

Catalog No: #AP82334



Package Size: #AP82334-1 50ug #AP82334-2 100ug #AP82334-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chymotrypsinogen 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	CTRB
Accession No.	Uniprot ID:P04813
Uniprot	P04813
GeneID	479649;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CTRB1

Note: This product is for in vitro research use only