

C-X-C motif chemokine 10

Catalog No: #AP82335



Package Size: #AP82335-1 50ug #AP82335-2 100ug #AP82335-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Small-inducible cytokine B10
Accession No.	Uniprot ID:Q5KSV9
Uniprot	Q5KSV9
GeneID	478432;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CXCL10

Note: This product is for in vitro research use only