

# Protein CYR61

Catalog No: #AP82336



Package Size: #AP82336-1 50ug #AP82336-2 100ug #AP82336-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein CYR61
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CCN family member 1,Cysteine-rich angiogenic inducer 61,Insulin-like growth factor-binding protein 10,IBP-10,Ccn1,Igfbp10
Accession No.	Uniprot ID:Q9ES72
Uniprot	Q9ES72
GeneID	83476;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cyr61

Note: This product is for in vitro research use only