

TLK1 (Phospho-Ser764) Antibody

Catalog No: #12419

Package Size: #12419-1 50ul #12419-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TLK1 (Phospho-Ser764) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms
Specificity	TLK1 (Phospho-Ser764) Antibody detects endogenous levels of TLK1 only when phosphorylated at Ser764
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TLK1 (Phospho-Ser764)
Target Name	TLK1
Modification	Phospho
Other Names	TLK1 (Phospho-Ser764) Antibody
Accession No.	Swiss-Prot#: Q9Y4F6NCBI Gene ID: 9874
Uniprot	Q9UKI8
GeneID	9874;
Target Species	human
Calculated MW	86kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only