

CX3C chemokine receptor 1

Catalog No: #AP82360



Package Size: #AP82360-1 50ug #AP82360-2 100ug #AP82360-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CX3C chemokine receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C-X3-C CKR-1,Fractalkine receptor,RBS11,Rbs11
Accession No.	Uniprot ID:P35411
Uniprot	P35411
GeneID	171056;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cx3cr1

Note: This product is for in vitro research use only