

Dihydrofolate reductase

Catalog No: #AP82365



Package Size: #AP82365-1 50ug #AP82365-2 100ug #AP82365-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dihydrofolate reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:P00376
Uniprot	P00376
GeneID	508809;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DHFR

Note: This product is for in vitro research use only