

Carnitine O-palmitoyltransferase 1, liver isoform

Catalog No: #AP82371



Package Size: #AP82371-1 50ug #AP82371-2 100ug #AP82371-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carnitine O-palmitoyltransferase 1, liver isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CPT1-L,Carnitine O-palmitoyltransferase I, liver isoform,CPT I,Carnitine palmitoyltransferase 1A,Cpt-1,Cpt1
Accession No.	Uniprot ID:P97742
Uniprot	P97742
GeneID	12894;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cpt1a

Note: This product is for in vitro research use only