

Coenzyme Q-binding protein COQ10 homolog A, mitochondrial

Catalog No: #AP82372

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP82372-1 50ug #AP82372-2 100ug #AP82372-3 1mg

Description

Product Name	Coenzyme Q-binding protein COQ10 homolog A, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	UNQ6192/PRO20219
Accession No.	Uniprot ID:Q96MF6
Uniprot	Q96MF6
GeneID	93058;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COQ10A

Note: This product is for in vitro research use only