

C-X-C chemokine receptor type 1

Catalog No: #AP82375



Package Size: #AP82375-1 50ug #AP82375-2 100ug #AP82375-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	C-X-C chemokine receptor type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CXC-R1,High affinity interleukin-8 receptor A,IL-8R A,CD181,Il8ra
Accession No.	Uniprot ID:P70612
Uniprot	P70612
GeneID	54258;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cxcr1

Note: This product is for in vitro research use only