

# Tumor necrosis factor receptor superfamily member EDAR



Catalog No: #AP82381

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP82381-1 50ug #AP82381-2 100ug #AP82381-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor receptor superfamily member EDAR
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Anhidrotic ectodysplasin receptor 1,Downless,Ectodermal dysplasia receptor,Ectodysplasin-A receptor,DI
Accession No.	Uniprot ID:Q9R187
Uniprot	Q9R187
GeneID	13608;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Edar

Note: This product is for in vitro research use only