

EF-hand domain-containing protein D2

Catalog No: #AP82385



Package Size: #AP82385-1 50ug #AP82385-2 100ug #AP82385-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	EF-hand domain-containing protein D2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Swiprosin-1,Sws1
Accession No.	Uniprot ID:Q4FZY0
Uniprot	Q4FZY0
GeneID	298609;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Efh2

Note: This product is for in vitro research use only