

CD130/gp130 (Phospho-Ser782) Antibody

Catalog No: #12423

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD130/gp130 (Phospho-Ser782) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu Ms
Specificity	CD130/gp130 (Phospho-Ser782) Antibody detects endogenous levels of CD130/gp130 only when phosphorylated at Ser782
Immunogen Type	Peptide
Immunogen Description	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S782 of human CD130 protein.
Target Name	CD130/gp130
Modification	Phospho
Other Names	IL6ST, CD130, CD130 antigen, IL-6R subunit beta, IL-6 receptor subunit beta, Membrane glycoprotein gp130, Membrane glycoprotein 130, IL-6RB, CDW130, gp130, IL-6R-beta
Accession No.	Swiss-Prot#: P40189NCBI Gene ID: 3572
Uniprot	P40189
GeneID	3572;
Target Species	human
Calculated MW	130kDa
Concentration	1.0mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of CD130 (pS782) expression in HEK293T (A), SGC7901 (B), LO2 (C), PC12 (D), AML12 (E) whole cell lysates.

Note: This product is for in vitro research use only