

Glutamate decarboxylase 1

Catalog No: #AP82392



Package Size: #AP82392-1 50ug #AP82392-2 100ug #AP82392-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamate decarboxylase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	67 kDa glutamic acid decarboxylase, GAD-67, Glutamate decarboxylase 67 kDa isoform, GAD67
Accession No.	Uniprot ID:A0PA85
Uniprot	A0PA85
GeneID	478794;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GAD1

Note: This product is for in vitro research use only