

Collagen alpha-1(VII) chain

Catalog No: #AP82396



Package Size: #AP82396-1 50ug #AP82396-2 100ug #AP82396-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen alpha-1(VII) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Long-chain collagen,LC collagen
Accession No.	Uniprot ID:Q63870
Uniprot	Q63870
GeneID	12836;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Col7a1

Note: This product is for in vitro research use only