

Desmoglein-1-alpha

Catalog No: #AP82398



Package Size: #AP82398-1 50ug #AP82398-2 100ug #AP82398-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Desmoglein-1-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Desmoglein-1,DG1,DGI,Desmosomal glycoprotein I,Dsg1
Accession No.	Uniprot ID:Q61495
Uniprot	Q61495
GeneID	13510;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dsg1a

Note: This product is for in vitro research use only