

C-reactive protein

Catalog No: #AP82403



Package Size: #AP82403-1 50ug #AP82403-2 100ug #AP82403-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	C-reactive protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	PTX1
Accession No.	Uniprot ID:O19062
Uniprot	O19062
GeneID	396842;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CRP

Note: This product is for in vitro research use only