

Elongation factor 1-alpha 1

Catalog No: #AP82406

Package Size: #AP82406-1 50ug #AP82406-2 100ug #AP82406-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Elongation factor 1-alpha 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EF-1-alpha-1, Elongation factor Tu, EF-Tu, Eukaryotic elongation factor 1 A-1, eEF1A-1, Leukocyte receptor cluster member 7, EEF1A, EF1A, LENG7
Accession No.	Uniprot ID:P68104
Uniprot	P68104
GeneID	1915;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: EEF1A1

Note: This product is for in vitro research use only