

C-X-C chemokine receptor type 5

Catalog No: #AP82408



Package Size: #AP82408-1 50ug #AP82408-2 100ug #AP82408-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C chemokine receptor type 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CXC-R5,Burkitt lymphoma receptor 1,Monocyte-derived receptor 15,MDR-15,CD185,BLR1,MDR15
Accession No.	Uniprot ID:P32302
Uniprot	P32302
GeneID	643;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CXCR5

Note: This product is for in vitro research use only