

C-X-C motif chemokine 14

Catalog No: #AP82409



Package Size: #AP82409-1 50ug #AP82409-2 100ug #AP82409-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	B-cell and monocyte-activating chemokine, Chemokine BRAK, Kidney-expressed chemokine CXC, MIP-2G, Small-inducible cytokine B14, Bmac, Kec, Ks1, Mip2g, Scyb14
Accession No.	Uniprot ID: Q9WUQ5
Uniprot	Q9WUQ5
GeneID	57266;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cxcl14

Note: This product is for in vitro research use only