

Cystatin-S

Catalog No: #AP82410



Package Size: #AP82410-1 50ug #AP82410-2 100ug #AP82410-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cystatin-S
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cystatin-1,Protein LM,Cyss
Accession No.	Uniprot ID:P19313
Uniprot	P19313
GeneID	296234;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cst4

Note: This product is for in vitro research use only