

ATP-dependent RNA helicase DDX42

Catalog No: #AP82412



Package Size: #AP82412-1 50ug #AP82412-2 100ug #AP82412-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP-dependent RNA helicase DDX42
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DEAD box protein 42
Accession No.	Uniprot ID:Q810A7
Uniprot	Q810A7
GeneID	72047;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ddx42

Note: This product is for in vitro research use only