

C-X-C chemokine receptor type 7

Catalog No: #AP82415



Package Size: #AP82415-1 50ug #AP82415-2 100ug #AP82415-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C chemokine receptor type 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CXC-R7,CXCR-7,Chemokine orphan receptor 1,G-protein coupled receptor RDC1 homolog,RDC-1,Cmkr1,Rdc1
Accession No.	Uniprot ID:O89039
Uniprot	O89039
GeneID	84348;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cxcr7

Note: This product is for in vitro research use only