

Death domain-associated protein 6

Catalog No: #AP82418



Package Size: #AP82418-1 50ug #AP82418-2 100ug #AP82418-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Death domain-associated protein 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Daxx,hDaxx,ETS1-associated protein 1,EAP1,Fas death domain-associated protein,BING2,DAP6
Accession No.	Uniprot ID:Q9UER7
Uniprot	Q9UER7
GeneID	1616;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DAXX

Note: This product is for in vitro research use only