

# Death-inducer obliterator 1

Catalog No: #AP82419

Package Size: #AP82419-1 50ug #AP82419-2 100ug #AP82419-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Death-inducer obliterator 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DIO-1,Death-associated transcription factor 1,DATF-1,Datf1,Dio1
Accession No.	Uniprot ID:Q8C9B9
Uniprot	Q8C9B9
GeneID	23856;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Dido1

Note: This product is for in vitro research use only