

Elongation factor 1-alpha 1

Catalog No: #AP82423

Package Size: #AP82423-1 50ug #AP82423-2 100ug #AP82423-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Elongation factor 1-alpha 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	EF-1-alpha-1, Elongation factor Tu, EF-Tu, Eukaryotic elongation factor 1 A-1, eEF1A-1, Eef1a
Accession No.	Uniprot ID:P10126
Uniprot	P10126
GeneID	13627;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Eef1a1

Note: This product is for in vitro research use only