

CX3C chemokine receptor 1

Catalog No: #AP82429



Package Size: #AP82429-1 50ug #AP82429-2 100ug #AP82429-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CX3C chemokine receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C-X3-C CKR-1,Beta chemokine receptor-like 1,CMK-BRL-1,CMK-BRL1,Fractalkine receptor,G-protein coupled receptor 13,V28,CMKBRL1,GPR13
Accession No.	Uniprot ID:P49238
Uniprot	P49238
GeneID	1524;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CX3CR1

Note: This product is for in vitro research use only