

Dipeptidyl peptidase 9

Catalog No: #AP82434



Package Size: #AP82434-1 50ug #AP82434-2 100ug #AP82434-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dipeptidyl peptidase 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DP9,Dipeptidyl peptidase IV-related protein 2,DPRP-2,Dipeptidyl peptidase IX,DPP IX,Dipeptidyl peptidase-like protein 9,DPLP9,DPRP2
Accession No.	Uniprot ID:Q86T12
Uniprot	Q86T12
GeneID	91039;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DPP9

Note: This product is for in vitro research use only