

Complement decay-accelerating factor

Catalog No: #AP82437



Package Size: #AP82437-1 50ug #AP82437-2 100ug #AP82437-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Complement decay-accelerating factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD55,CR,DAF
Accession No.	Uniprot ID:P08174
Uniprot	P08174
GeneID	1604;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD55

Note: This product is for in vitro research use only