

## Discoidin domain-containing receptor 2

Catalog No: #AP82438



Package Size: #AP82438-1 50ug #AP82438-2 100ug #AP82438-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Discoidin domain-containing receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Discoidin domain receptor 2,CD167 antigen-like family member B,Neurotrophic tyrosine kinase, receptor-related 3,Receptor protein-tyrosine kinase TKT,Tyrosine-protein kinase TYRO10,CD167b,Ntrkr3,Tkt,Tyro10
Accession No.	Uniprot ID:Q62371
Uniprot	Q62371
GeneID	18214;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ddr2

Note: This product is for in vitro research use only