

Beta-Ala-His dipeptidase

Catalog No: #AP82440

Package Size: #AP82440-1 50ug #AP82440-2 100ug #AP82440-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-Ala-His dipeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CNDP dipeptidase 1,Carnosine dipeptidase 1,Cn1
Accession No.	Uniprot ID:Q66HG3
Uniprot	Q66HG3
GeneID	307212;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cndp1

Note: This product is for in vitro research use only