

Cytosolic non-specific dipeptidase

Catalog No: #AP82441



Package Size: #AP82441-1 50ug #AP82441-2 100ug #AP82441-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytosolic non-specific dipeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CNDP dipeptidase 2,Carnosine dipeptidase II,Epididymis secretory protein Li 13,Glutamate carboxypeptidase-like protein 1,Peptidase A,CN2,CPGL,HEL-S-13,PEPA
Accession No.	Uniprot ID:Q96KP4
Uniprot	Q96KP4
GeneID	55748;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CNDP2

Note: This product is for in vitro research use only