

D(1A) dopamine receptor

Catalog No: #AP82442

Package Size: #AP82442-1 50ug #AP82442-2 100ug #AP82442-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	D(1A) dopamine receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dopamine D1 receptor,Drd1a,Gpcr15
Accession No.	Uniprot ID:Q61616
Uniprot	Q61616
GeneID	13488;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Drd1

Note: This product is for in vitro research use only