

E3 ubiquitin-protein ligase DTX1

Catalog No: #AP82443



Package Size: #AP82443-1 50ug #AP82443-2 100ug #AP82443-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase DTX1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FXI-T1, Protein deltex-1, Deltex1, RING-type E3 ubiquitin transferase DTX1
Accession No.	Uniprot ID:Q61010
Uniprot	Q61010
GeneID	14357;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Dtx1

Note: This product is for in vitro research use only