

Deoxyuridine 5'-triphosphate nucleotidohydrolase, mitochondrial

Catalog No: #AP82444

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP82444-1 50ug #AP82444-2 100ug #AP82444-3 1mg

Description

Product Name	Deoxyuridine 5'-triphosphate nucleotidohydrolase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	dUTPase,dUTP pyrophosphatase
Accession No.	Uniprot ID:P33316
Uniprot	P33316
GeneID	1854;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DUT

Note: This product is for in vitro research use only