

Dickkopf-like protein 1

Catalog No: #AP82445

Package Size: #AP82445-1 50ug #AP82445-2 100ug #AP82445-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dickkopf-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cancer/testis antigen 34,CT34,Protein soggy-1,SGY-1,SGY1,UNQ735/PRO1429
Accession No.	Uniprot ID:Q9UK85
Uniprot	Q9UK85
GeneID	27120;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DKKL1

Note: This product is for in vitro research use only