

ATP-dependent RNA helicase DDX42

Catalog No: #AP82447



Package Size: #AP82447-1 50ug #AP82447-2 100ug #AP82447-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP-dependent RNA helicase DDX42
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DEAD box protein 42, RNA helicase-like protein, RHELP, RNA helicase-related protein, RNAHP, SF3b DEAD box protein, Splicing factor 3B-associated 125 kDa protein, SF3b125
Accession No.	Uniprot ID:Q86XP3
Uniprot	Q86XP3
GeneID	11325;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DDX42

Note: This product is for in vitro research use only