

C-X-C motif chemokine 3

Catalog No: #AP82451



Package Size: #AP82451-1 50ug #AP82451-2 100ug #AP82451-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytokine-induced neutrophil chemoattractant 2,CINC-2,Macrophage inflammatory protein 2-alpha/beta,MIP2-alpha/beta,Cinc2
Accession No.	Uniprot ID:Q10746
Uniprot	Q10746
GeneID	171551;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cxcl3

Note: This product is for in vitro research use only