

Dipeptidyl peptidase 4

Catalog No: #AP82452



Package Size: #AP82452-1 50ug #AP82452-2 100ug #AP82452-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dipeptidyl peptidase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Activation molecule 3,ACT3,Adenosine deaminase complexing protein,ADCP-I,Dipeptidyl peptidase IV,DPP IV,T-cell activation antigen CD26,WC10,CD26,CD26
Accession No.	Uniprot ID:P81425
Uniprot	P81425
GeneID	281122;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DPP4

Note: This product is for in vitro research use only